

Mag-Gro™

MAGNETIC HAIR AND SCALP STIMULATOR

Patent Pending

It's All Natural

It's Easy to Understand

And it works - GUARANTEED

Every product that claims to slow or stop hair loss does so by attempting to do only one thing:

Increase blood flow to the scalp.

- Minoxidil based products – vaso dilation
- Vibrating hair brush – scalp massage
- Laser comb – skin stimulation

Mag-Gro increases blood flow to the scalp by using a three-step process:

- Compression of skin and hair follicles
- Physical scalp massage & stimulation
- Rotating magnetic field enhances blood flow

Mag-Gro can increase blood flow *better, faster and cheaper than any other product in use today.* But Mag-Gro may do more than stop hair loss.

Mag-Gro can also create a healthier scalp and better-looking hair *in three additional ways:*

1. Generate full size hair in bald areas
2. Increase hair oil for more softness and shine
3. Put color back into gray hair and prevent gray hair from forming

Mag-Gro uses an all-natural process to achieve these results: no chemicals on your scalp and no hormones in your body.

MAG-GRO MAKES SENSE

HOW DOES MAG-GRO WORK?

Mag-Gro is a massaging tool that stimulates blood flow in three ways:

1. Compression of the individual hair follicles by the specially designed little fingers;
2. Massaging action of the roller while the user presses down gently on the scalp
3. Four magnets are embedded in the roller. The poles are opposing which draws blood back and forth through the capillaries.

Mag-Gro may help your scalp's natural hair growth process in 4 ways:

- Increased blood flow through existing capillaries and/or
- Stimulation of new capillary growth
- Increased circulation and/or distribution of hair oil from the sebaceous gland
- Impact on the hair's melanocyte cells which increases melanin into the hair shaft. This will decrease or eliminate existing gray hair and stop new gray hair from forming.



Mag-Gro should be used on a regular basis in order to maintain a healthy head of hair. Just as *you brush your teeth to make them look better and last longer, so should you use Mag-Gro to make your scalp and hair look and feel better.*

DIRECTIONS: Insert one or two fingers in each end of the roller. Place the roller on your scalp. Gently guide the Roller back and forth over the scalp.

The back and forth motion will cause the specially designed fingers to press down on the scalp and the individual hair follicles. This tends to stimulate activity in the hair follicle. You should begin to feel blood flowing and a general tingling sensation in the scalp.

Press down gently as you roll Mag-Gro over the scalp. Do not allow the roller to skid across the skin. In time you will be able to press down more firmly.

Roll Mag-Gro over the entire scalp including the back of the head, the sides of the head, the top, temples and even the forehead. Stimulate the entire scalp.

The opposing magnets create a magnetic field around Mag-Gro. This rotating magnetic field will further enhance blood flow to the hair follicle.

After rolling your scalp for 30-45 seconds the skin will appear pink and a tingling sensation can be felt. This is an indication that blood is returning to your scalp.

Try to use Mag-Gro at least once a day. Use it in the shower. Use it during a commercial on TV. Keep one Mag-Gro in your desk at work, or in your gym bag.

Use Mag-Gro as often as you like but let your scalp's sensitivity be your guide. The first time user should be careful not to over stimulate the scalp or abrade the skin. At the same time, the more blood that gets to the follicle, the better it is for your scalp and hair.

Make Mag-Gro part of your normal health maintenance routine:

- Wash your face,
- Brush your teeth, and

- Roll your scalp with Mag-Gro

For a more detailed explanation please go to www.rollyourscalp.com



WARNING: Move Mag-Gro slowly and deliberately across the scalp. DO NOT use a rapid motion. Press down on the roller only to tolerance. Allow Mag-Gro to roll across the scalp. DO NOT cause the roller to scrap or skid across the skin.

It is not possible to damage the scalp or cause discomfort unless you press down too hard on the roller. If you feel discomfort at any time, STOP what you are doing immediately.

Wait a minute or two before continuing to use Mag-Gro. When you continue, press down on the scalp with less force.

DO NOT have unrealistic expectations. Mag-Gro is an all-natural process. If you are bald, Mag-Gro will not give you a full head of hair overnight. Mag-Gro will, however, do more in the long run to improve the condition of your hair and scalp than any other product on the market today.

GUARANTEE: If Mag-Gro is used as directed at least once daily for 90 days, you are guaranteed to experience a healthier scalp and better looking hair or you can request a full refund by calling,

1 877 ROLL SCALP (765-5722)